




Meeka Microscope Safety

*Congratulations, you now have microscopes at your school! You and your students are going to be amazed at what you will see.

Microscope Cart Use Safety Rules:

To be safe when using the microscope cart you must follow these safety rules to avoid injury or damage due to loss of control.

1. **Use two adults or students when wheeling the cart up or down inclined surfaces** such as a ramp to your classroom. Have one student or adult on the side of the cart and the other on the downhill side to support the cart as you descend a ramp. Do not try to take the cart on a ramp alone.
2. **Set the wheel brakes when the cart is not in use**, especially if the cart is near sloping ground. Wheel brakes are located on the two straight wheels and are labeled with brake on and brake off pedal positions.
3. Be aware that magnetic doors pull off of the cart front and are put on the side of the cart.



Enjoy using your
microscope cart!
Sincerely, Meeka
Microscope

